

medi

New in this guide:
mediven harmony – short

medi formulary guide FP10 – GP10

September 2016



medi. I feel better.



medi

medi is the largest manufacturer of compression garments in the world and has built a reputation for quality and innovation over the past 60 years. The main focus has always been on clinical benefits, comfort and concordance.

Improvements in materials, manufacturing techniques and knitting technology have been adopted by medi to help produce compression garments of the highest standard.

The end user is always the primary focus and constant feedback has brought continued innovation and product development.

All medi garments featured in this guide are available on the Drug Tariff in the UK and have been fully approved by the Department of Health as being suitable for patient prescribing. There are cost savings associated with the correct and timely use of compression garments, particularly when treating those with venous disease, lymphoedema and tissue viability issues associated with venous leg ulcers. This compact formulary guide is designed to assist with the correct selection of mediven and juxtafit garments.



Within medi UK, we are committed to supporting you in continuing your professional development and therefore have developed learning packages that have CPD accreditation.

Obtaining CPD accreditation means that our training has been assessed and scrutinised by the Professional CPD Association to ensure it is of the highest quality and integrity. medi UK recognises the importance of CPD and how this can benefit both your personal development and ultimately the care that your patients receive. medi UK want to ensure that we offer you the best training possible in order to help support you and your patients.

For information on our CPD accredited events – please visit <http://www.medi.biz/accredited-courses> or email CPD@mediuk.co.uk





medi World of Compression

A head start in treatment. Convincing healthcare concepts. Effective products.

medi World of Compression: signifies over 65 years competence in the field of compression. medi aim to promote the highest care and concordance to treatment – always focussing on the best interests of the patient. Together we can help to shape a future of successful treatment for you and your patients.

For clinicians:

- Nationwide clinical support, tailored in-house training and CPD courses
- Quality RAL and British Standard compression
- medi Therapy Chains for wound healing and chronic oedema / lymphoedema

We offer your patients:

- The complete care pathway from skin care & debridement to wound healing & confidence in managing their conditions
- Increased mobility
- More zest in life and a better quality of life

Compression standards

There are three compression standards commonly used in the UK, and practitioners need to be aware of the differences as detailed in the table below.

All mediven garments featured in this guide are knitted to the highest quality RAL standard. Certification at this level requires stringent testing by the Hohenstein Institute and is the gold standard for compression garments. The European standard was considered but then discarded in 2002 as there was no testing protocol in place.

Indications for compression garments – effective quality tested compression:

RAL class 1 (18 – 21mmHg)

Prevention or early treatment of lymphoedema, venous disease and for use on palliative limbs.

RAL class 2 (23 – 32mmHg)

Mild/moderate lymphatic conditions, pronounced varicosis with oedema present and after healing minor ulcers.

RAL class 3 (34 – 46mmHg)

Lymphoedema, significant varicosis or after healing long standing ulcers; use if condition fails to respond to class 2.

RAL Standard Hohenstein (HOSY) tested

Contra-indications:

According to the guidelines of the German Phlebological Society (published in 1996, German Journal “Phlebologie”).

Absolute contra-indications:

Advanced, peripheral arterial occlusive disease, decompensated heart disease, septic phlebitis, phlegmasia coerulea dolens.

Relative contra-indications:

Suppurative dermatoses, intolerance of compression stocking fabric, sensory disturbances of the limb, advanced peripheral neuropathy, primary chronic arthritis. medi assumes no liability as a result of any contra-indicated use of this product.

mediven® offers many unique features

Optimising wearing comfort



Clima Comfort

Clima Comfort ensures maximum comfort

- Cool in summer, warm in winter
- Moisture wicking
- Breathable



Soft Elastic

Soft Elastic

- Improved wearing comfort due to the noticeably increased toe space for closed toe options



Clima Fresh

Clima Fresh for more freshness

- Antibacterial properties built into the knit
- Odour control
- Unique spinning technique that guarantees the effect, even after daily washing



Perfect Fit

Perfect Fit

- Very precise compression
- Wide variety of off-the-shelf sizes



Double ellipses

Double Ellipses, innovative fabric

- The double ellipses on the point and in the crease of the elbow offer comfort, and freedom of movement
- Innovative soft fabric – offering effective compression



Variety

Made-to-measure options for non-standard sizes

- For severe indications
- For any leg shapes
- Special features / options available



Oeko-Tex Standard 100

All compression stockings and garments, whether circular-knit or flat-knit, thrombosis prophylaxis stockings and supports are regularly examined and rated according to strict standards by independent test centres – eg. Institute Hohenstein – for any toxic substances potentially damaging to human health. Granting of the Oeko-Tex-Standard 100 certificate guarantees: tested quality for the protection of your health.



Very good skin tolerability and suitability for everyday wear

- Without cotton or rubber

mediven product features

Product features	mediven plus	mediven elegance	mediven for men	mediven active	mediven mondi	mediven ulcer kit	mediven esprit	mediven harmony	duomed soft
 Clima Comfort	●	●	●	●	●		●	●	
 Clima Fresh	●	●	●		●		●		
 Soft Elastic		●							
 Perfect Fit	●	●	●	●	●		●	●	

Product features	mediven plus	mediven elegance	mediven for men	mediven active	mediven mondi	mediven ulcer kit	mediven esprit	mediven harmony	duomed soft
 Variety	●	●			●		●	●	
 Double ellipses								●	
 Anatomic Fit					●		●		
 LATEX	●	●	●	●	●		●		●

mediven plus®

The open toe RAL compression stocking, in classes 1, 2 and 3.

Features

- the plus-heel with its Y-shape provides for easy donning, a secure fit and a free feeling in the instep area - all day long
- the 5 cm wide cuff on the knee high socks ensures maximum comfort and a secure hold
- accurate graduated compression
- wide range of variants available – less made-to-measure required
- Clima Comfort technology ensures rapid evaporation of moisture



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

Product benefits



Suitable for use with this product and available on prescription: **medi 2in1**

helps to put on and take off open toe and closed toe compression garments. Product code DT904, PIP code 329-8932



mediven elegance®

The fashionable stocking with a lace topband.

Features

- a noticeably lower sense of pressure with full medical compression efficacy
- a perfect fit and outstanding comfort in wear
- fashionable elegance and beauty
- the below-knee stocking cuff is particularly wide and soft, which gives a perfect hold and makes the stocking so comfortable to wear



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

Product benefits



Suitable for use with this product and available on prescription:
medi 2in1

helps to put on and take off open toe and closed toe compression garments.
Product code DT904, PIP code 329-8932



mediven® for men

The elegant below knee sock for men.

Features

- 'Business Rib' design
- anatomically shaped toes
- Clima Comfort technology ensures rapid evaporation of moisture
- cool in summer, warm in winter
- machine wash and tumble dry on cool settings
- stylish choice of colours



Product benefits



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

Suitable for use with this product and available on prescription: medi 2in1

helps to put on and take off open toe and closed toe compression garments.
Product code DT904, PIP code 329-8932



mediven active®

The robust and hard-wearing below knee sock.

Features

- fine ribbed design
- extra wide cuff
- seamless toe section
- reinforced heel and comfort reinforced sole
- Clima Comfort technology ensures rapid evaporation of moisture
- cool in summer, warm in winter
- machine washable and tumble dry on cool settings
- 5 colours, 7 sizes and 2 lengths available. Petite length only available in black



Product benefits



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

Suitable for use with this product and available on prescription: medi 2in1

helps to put on and take off open toe and closed toe compression garments. Product code DT904, PIP code 329-8932



Benefits of using RAL compression hosiery, cost-effective patient-centred care.

Advances in technology over the last ten years have allowed for greater choice of compression therapy with proven clinical results and increased concordance.

A recent study by Caroline Dowsett, Nurse Consultant Tissue Viability, has shown that using RAL garments has benefits when treating and preventing recurrence of venous leg ulcers.* Healing rates were shown to improve from 36 % to 72 % at 12 weeks and from 40 % to 100 % at 24 weeks. Recurrence rates were reduced from 18-20 % to 5.8 %. This service was re-audited one year later and the recurrence rates had dropped further to 2.8 %.

The reduction in healing times and reduced recurrence rates has greatly reduced the cost of venous ulcer care for Newham PCT.

**Available free of charge from medi UK are:
copies of both papers quoted, a CD of Caroline's presentation or your copy of 'A Leg Ulcer Service Development toolkit' – to help you achieve and monitor cost savings and clinical effectiveness.
Please contact our customer service team on 01432 373500.**



*Dowsett, C. (2011): Treatment and prevention of recurrence of venous leg ulcers using RAL hosiery, in: Wounds UK, Vol. 7, no 1, p. 115-119
Dowsett, C (2012): Modernising leg ulcer services through preventing recurrence. Wounds UK, Vol 8, No.1

Useful tips for measuring:

Below Knee garments:

- 1) measure the circumference at the narrowest part of the ankle ('B'). If there is oedema present, and there is no obvious narrow part, use the point 2 fingers width up from the tip of the lateral malleolus
 - 2) measure the circumference at the widest part of the calf ('C')
 - 3) measure the length of the limb at the back from the heel (flat on the floor) to a point 2 finger widths down from the knee flexure
- 1 and 2 above give you the measurements to decide the size of the garment (e.g.3 or 4). 3 helps you to select the correct length of garment (standard or petite).

For example:

circumference B=24 cm, circumference C=40 cm, therefore could be a size 3 or 4:

Circumference: size in cm				
Size	II	III	IV	V
⌂C	30 – 37	33 – 40	35 – 43	37 – 46
⌂B	20 – 22	22 – 24	24 – 26	26 – 28

If the measurements fall between two sizes, use your clinical judgement for each patient. The larger size would be appropriate for an older person who will struggle with compression, or the smaller size for someone who is young and/or active.

Thigh Length garments:


Measure 'B' and 'C' as for Below Knee garments and take an additional circumference measurement at 'G' (thigh at point where top of garment will be). The length is from the heel to 4 fingers width down from the inguinal crease, this determines whether the patient requires a standard or petite length garment.

measuring for mediven mondi-made-to-measure.

Please contact customer service team on 01432 373500 and ask for 'FP10 m2m brochure'.

Measurement guide leg

Length: size in cm	Petite	Regular
AD calf stocking	34 – 38	39 – 44
mediven active only AD calf sock	37 – 42	43 – 49
AG thigh stocking	62 – 71	72 – 83



Circumference: size in cm

Size	I	II	III	IV	V	VI	VII
⌘G thigh with silicone topband extra wide	49 – 57	53 – 62	57 – 67	61 – 72	65 – 77	69 – 81	73 – 85
⌘G thigh with silicone topband	43 – 48	45 – 52	49 – 56	53 – 60	56 – 64	60 – 68	64 – 72
⌘G thigh with waist attachment, tights	43 – 57	45 – 62	49 – 67	53 – 72	56 – 77	60 – 81	64 – 85
⌘C calf extra wide	34 – 38	37 – 41	40 – 44	42 – 46	44 – 48	46 – 50	48 – 52
⌘C calf	28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
⌘B ankle	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32

RAL compression class 1

Suitable for treatment of mild lymphoedema, mild varicosis without oedema and for early varicosis in pregnancy.

Product codes for compression class 1 leg garments / 18 – 21mmHg

Styles	open toe	closed toe	AD	AD pet.	AD	AD pet.	AG /SN	AG pet./SN	AG /SN	AG pet./SN
										
			calf	calf petite	calf x-wide	calf x-wide petite	thigh with topband	thigh petite with topband	thigh x-wide with topband	thigh x-wide petite with topband
mediven plus	•		DT101	DT102	DT103	DT104	DT108	DT106	DT108X	DT106X
mediven elegance		•	DT185	DT186			DT189	DT190	DT189X	DT190X
mediven for men		•	DT342	DT343						
mediven active		•	DT180	DT181						

Product codes for compression class 1 leg garments / 18 – 21mmHg

Styles	open toe	closed toe	AG	AG pet.	AT	AT pet.	beige	black	navy	grey	white	brown
												
			thigh with waist attachment	thigh petite with waist attachment	tights	tights petite						
			right	left	right	left						
mediven plus	•		DT117	DT119	DT118	DT120	•	•				
mediven elegance		•			DT191	DT192	•	•				
mediven for men		•						•	•	•		
mediven active		•						•	•	•	•	•

mediven plus® PIP codes ccl 1

ccl 1 mediven plus PIP codes below knee									
		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT101/0/...	215-6016	323-5405	323-5413	323-5421	323-5439	323-5447	323-5454
	black	DT101/5/...	323-5538	323-5546	323-5553	323-5561	323-5579	323-5587	323-5595
open toe, petite	beige	DT102/0/...	323-5678	323-5686	323-5694	323-5702	323-5710	323-5728	323-5736
	black	DT102/5/...	323-5819	323-5827	323-5835	323-5843	323-5850	323-5868	323-5876
open toe, extra wide	beige	DT103/0/...	323-5462	323-5470	323-5488	323-5496	323-5504	323-5512	323-5520
	black	DT103/5/...	323-5603	323-5611	323-5629	323-5637	323-5645	323-5652	323-5660
open toe, petite, extra wide	beige	DT104/0/...	323-5744	323-5751	323-5769	323-5777	323-5785	323-5793	323-5801
	black	DT104/5/...	323-5884	323-5892	323-5900	323-5918	323-5926	323-5934	323-5942

ccl 1 mediven plus PIP codes thigh length

Silicone topband		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT108/0/...	323-7492	323-7500	323-7518	323-7526	323-7534	323-7542	323-7559
	black	DT108/5/...	323-7567	323-7575	323-7583	323-7591	323-7609	323-7617	323-7625
open toe, thigh extra wide	beige	DT108X/0/...	334-0429	334-0437	334-0445	334-0452	334-0460	334-0478	334-0486
open toe, petite	beige	DT106/0/...	323-6916	323-6924	323-6932	323-6940	323-6957	323-6965	323-6973
	black	DT106/5/...	323-6981	323-6999	323-7005	323-7013	323-7021	323-7039	323-7047
open toe, thigh extra wide, petite	beige	DT106X/0/...	334-0569	334-0577	334-0585	334-0593	334-0601	334-0619	334-0643
Waist attachment		medi code	I	II	III	IV	V	VI	VII
open toe, right	beige	DT117/0/...	323-7815	323-7823	323-7831	323-7849	323-7856	323-7864	323-8011
	black	DT117/5/...	323-8029	323-8037	323-8045	323-8052	323-8060	323-8078	323-8086
open toe, petite, right	beige	DT118/0/...	323-7195	323-7203	323-7211	323-7229	323-7237	323-7245	323-7252
	black	DT118/5/...	323-7260	323-7278	323-7286	323-7294	323-7302	323-7310	323-7328
open toe, left	beige	DT119/0/...	323-7633	323-7641	323-7658	323-7666	323-7674	323-7682	323-7690
	black	DT119/5/...	323-7708	323-7716	323-7724	323-7732	323-7740	323-7757	323-7765
open toe, petite, left	beige	DT120/0/...	323-7054	323-7062	323-7070	323-7088	323-7096	323-7104	323-7112
	black	DT120/5/...	323-7120	323-7138	323-7146	323-7153	323-7161	323-7179	323-7187

mediven elegance® PIP codes ccl 1

ccl 1 mediven elegance PIP codes below knee									
		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT185/0/...	338-8964	338-8972	338-8980	338-8998	338-9004	338-9012	338-9020
	black	DT185/5/...	338-9095	338-9087	338-9079	338-9061	338-9053	338-9046	338-9038
closed toe, petite	beige	DT186/0/...	338-9103	338-9145	338-9152	338-9160	338-9178	338-9186	338-9194
	black	DT186/5/...	338-9202	338-9210	338-9228	338-9236	338-9244	338-9251	338-9269

ccl 1 mediven elegance PIP codes thigh length/tights

Silicone topband		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT189/0/...	338-9277	338-9285	338-9293	338-9301	338-9319	338-9327	338-9335
	black	DT189/5/...	338-9665	338-9673	338-9681	338-9699	338-9707	338-9715	338-9723
closed toe, thigh extra wide	beige	DT189X/0/...	338-9590	338-9608	338-9616	338-9624	338-9632	338-9640	338-9657
closed toe, petite	beige	DT190/0/...	338-9806	338-9814	338-9822	338-9830	338-9848	338-9855	338-9863
	black	DT190/5/...	338-9871	338-9889	338-9897	338-9905	338-9913	338-9921	338-9939
closed toe, thigh extra wide, petite	beige	DT190X/0/...	338-9947	338-9954	338-9962	338-9970	338-9988	338-9996	339-0002
Tights		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT191/0/...	333-8944	333-8951	333-8969	333-8977	333-8985	333-8993	333-9009
	black	DT191/5/...	333-9017	333-9025	333-9033	333-9041	333-9058	333-9066	333-9074
closed toe, petite	beige	DT192/0/...	333-9082	333-9090	333-9108	333-9116	333-9124	333-9132	333-9140
	black	DT192/5/...	333-9157	333-9165	333-9173	333-9181	333-9207	333-9215	333-9223

mediven® for men PIP codes ccl 1

ccl 1 mediven for men PIP codes									
		medi code	I	II	III	IV	V	VI	VII
closed toe, below knee	black	DT342/5/...	356-8433	349-6734	349-6742	349-6759	349-6767	349-6775	356-8441
	navy	DT342/3/...	356-8417	349-6684	349-6692	349-6700	349-6718	349-6726	356-8425
	grey	DT342/2/...	360-2885	360-2893	360-2901	360-2919	360-2927	360-2935	360-2943
closed toe, below knee, petite	black	DT343/5/...	356-8474	349-6833	349-6841	349-6858	349-6866	349-6874	356-8482
	navy	DT343/3/...	356-8458	349-6783	349-6791	349-6809	349-6817	349-6825	356-8466
	grey	DT343/2/...	360-2950	360-2968	360-2976	360-2984	360-2992	360-3008	360-3016

mediven active® PIP codes ccl 1

ccl 1 mediven active PIP codes									
		medi code	I	II	III	IV	V	VI	VII
standard	black	DT180/5/...	376-6466	376-6474	376-6482	376-6490	376-6508	376-6516	376-6524
	brown	DT180/6/...	376-6607	376-6615	376-6623	376-6631	376-6649	376-6656	376-6664
	grey	DT180/2/...	376-6326	376-6334	376-6342	376-6359	376-6367	376-6375	376-6383
	navy	DT180/3/...	376-6391	376-6409	376-6417	376-6425	376-6433	376-6441	376-6458
	white	DT180/7/...	367-6532	376-6540	376-6557	376-6565	376-6573	376-6581	376-6599
petite	black	DT181/5/...	367-6813	367-6821	376-6839	376-6847	367-6854	376-6862	376-6870

RAL compression class 2

Suitable for mild/moderate lymphatic conditions, pronounced varicosis with oedema present, after healing minor ulcers.

Product codes for compression class 2 leg garments / 23 – 32mmHg

Styles	open toe	closed toe	AD	AD pet.	AD	AD pet.	AG /SN	AG pet./SN	AG /SN	AG pet./SN
										
			calf	calf petite	calf x-wide	calf x-wide petite	thigh with topband	thigh petite with topband	thigh x-wide with topband	thigh x-wide petite with topband
mediven plus	•		DT201	DT202	DT203	DT 204	DT 200	DT199	DT200X	DT199X
mediven elegance		•	DT285	DT286			DT289	DT290	DT289X	DT290X
mediven for men		•	DT340	DT341						
mediven active		•	DT280	DT281						

Product codes for compression class 2 leg garments / 23 – 32mmHg

Styles	open toe	closed toe	AG		AG pet.		AT	AT pet.	beige	black	navy	grey	white	brown
														
			thigh with waist attachment		thigh petite with waist attachment		tights	tights petite						
			right	left	right	left								
mediven plus	•		DT207	DT209	DT208	DT210			•	•				
mediven elegance		•					DT291	DT292	•	•				
mediven for men		•								•	•	•		
mediven active		•								•	•	•	•	•

mediven plus® PIP codes ccl 2

ccl 2 mediven plus PIP codes below knee									
		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT201/0/...	033-5752	323-5959	323-5967	323-5975	323-5983	323-5991	323-6007
	black	DT201/5/...	323-6080	323-6098	323-6106	323-6114	323-6122	323-6130	323-6148
open toe, petite	beige	DT202/0/...	323-6239	323-6247	323-6254	323-6262	323-6270	323-6288	323-6296
	black	DT202/5/...	323-6379	323-6387	323-6395	323-6403	323-6411	323-6429	323-6437
open toe, extra wide	beige	DT203/0/...	215-6024	323-6023	323-6031	323-6049	323-6056	323-6064	323-6072
	black	DT203/5/...	323-6155	323-6163	323-6171	323-6189	323-6197	323-6213	323-6221
open toe, petite, extra wide	beige	DT204/0/...	323-6304	323-6312	323-6320	323-6338	323-6346	323-6353	323-6361
	black	DT204/5/...	323-6445	323-6452	323-6460	323-6478	323-6486	323-6494	323-6502

ccl 2 mediven plus PIP codes thigh length

Silicone topband		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT200/0/...	323-8797	323-8805	323-8813	323-8821	323-8839	323-8847	323-8854
	black	DT200/5/...	323-8862	323-8870	323-8888	323-8896	323-8904	323-8912	323-8920
open toe, thigh extra wide	beige	DT200X/0/...	334-1385	334-1393	334-1401	334-1419	334-1427	334-1435	334-1443
open toe, petite	beige	DT199/0/...	323-8235	323-8243	323-8250	323-8268	323-8276	323-8284	323-8292
	black	DT199/5/...	323-8300	323-8318	323-8326	323-8334	323-8342	323-8359	323-8367
open toe, thigh extra wide, petite	beige	DT199X/0/...	334-1658	334-1666	334-1674	334-1682	334-1690	334-1708	334-1716
Waist attachment		medi code	I	II	III	IV	V	VI	VII
open toe, right	beige	DT207/0/...	323-9076	323-9084	323-9092	323-9100	323-9118	323-9126	323-9134
	black	DT207/5/...	323-9142	323-9159	323-9167	323-9175	323-9183	323-9191	323-9209
open toe, petite, right	beige	DT208/0/...	323-8516	323-8524	323-8532	323-8540	323-8557	323-8565	323-8573
	black	DT208/5/...	323-8581	323-8599	323-8607	323-8615	323-8623	323-8631	323-8649
open toe, left	beige	DT209/0/...	323-8938	323-8946	323-8953	323-8961	323-8979	323-8987	323-8995
	black	DT209/5/...	323-9001	323-9019	323-9027	323-9035	323-9043	323-9050	323-9068
open toe, petite, left	beige	DT210/0/...	323-8375	323-8383	323-8391	323-8409	323-8417	323-8425	323-8433
	black	DT210/5/...	323-8441	323-8458	323-8466	323-8474	323-8482	323-8490	323-8508

mediven elegance® PIP codes ccl 2

ccl 2 mediven elegance PIP codes below knee									
		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT285/0/...	339-0408	339-0416	339-0424	339-0432	339-0440	339-0457	339-0465
	black	DT285/5/...	339-0481	339-0499	339-0507	339-0515	339-0523	339-0531	339-0549
closed toe, petite	beige	DT286/0/...	339-0572	339-0580	339-0598	339-0606	339-0614	339-0622	339-0648
	black	DT286/5/...	339-0655	339-0671	339-0689	339-0705	339-0713	339-0721	339-0739

ccl 2 mediven elegance PIP codes thigh length/tights

Silicone topband		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT289/0/...	339-0770	339-0796	339-0820	339-0838	339-0846	339-0853	339-0861
	black	DT289/5/...	339-0879	339-0887	339-0895	339-0903	339-0911	339-0929	339-0937
closed toe, thigh extra wide	beige	DT289X/0/...	339-1976	339-1984	339-1992	339-2008	339-2016	339-2024	339-2032
closed toe, petite	beige	DT290/0/...	339-2115	339-2123	339-2131	339-2149	339-2156	339-2164	339-2172
	black	DT290/5/...	339-2180	339-2198	339-2206	339-2214	339-2222	339-2230	339-2248
closed toe, thigh extra wide, petite	beige	DT290X/0/...	339-2255	339-2263	339-2271	339-2289	339-2297	339-2305	339-2313
Tights		medi code	I	II	III	IV	V	VI	VII
closed toe	beige	DT291/0/...	333-9330	333-9348	333-9355	333-9363	333-9371	333-9389	333-9397
	black	DT291/5/...	333-9405	333-9413	333-9421	333-9439	333-9447	333-9454	333-9462
closed toe, petite	beige	DT292/0/...	333-9470	333-9488	333-9496	333-9504	333-9512	333-9520	333-9538
	black	DT292/5/...	333-9546	333-9553	333-9561	333-9579	333-9587	333-9595	333-9603

mediven® for men PIP codes ccl 2

ccl 2 mediven for men PIP codes									
		medi code	I	II	III	IV	V	VI	VII
closed toe, below knee	black	DT340/5/...	360-4451	334-0791	334-0809	334-0817	334-0825	334-0833	360-4469
	navy	DT340/3/...	360-4436	334-0742	334-0759	334-0767	334-0775	334-0783	360-4444
	grey	DT340/2/...	360-3024	360-3032	360-3040	360-3057	360-3065	360-3073	360-3081
closed toe, below knee, petite	black	DT341/5/...	360-4493	334-0841	334-0858	334-0866	334-0874	334-0882	360-4501
	navy	DT341/3/...	360-4477	334-0890	334-0916	334-0932	334-0940	334-0957	360-4485
	grey	DT341/2/...	360-3099	360-3107	360-3115	360-3123	360-3131	360-3149	360-3156

mediven active® PIP codes ccl 2

ccl 2 mediven active PIP codes									
		medi code	I	II	III	IV	V	VI	VII
standard	black	DT280/5/...	376-7183	376-7209	376-7217	376-7225	376-7233	376-7241	376-7258
	brown	DT280/6/...	376-7332	376-7340	376-7357	376-7365	376-7373	376-7381	376-7399
	grey	DT280/2/...	376-7027	376-7035	376-7043	376-7050	376-7068	376-7076	376-7084
	navy	DT280/3/...	376-7092	376-7100	376-7118	376-7126	376-7134	376-7142	376-7167
	white	DT280/7/...	376-7266	376-7274	376-7282	376-7290	376-7308	376-7316	376-7324
petite	black	DT281/5/...	376-7829	376-7837	376-7845	376-7852	376-7860	376-7878	376-8025

Product codes for compression class 3 leg garments / 34 – 46mmHg

Styles	open toe	closed toe	AG		AG pet.		AT	AT pet.	beige	black
										
			thigh with waist attachment		thigh petite with waist attachment		tights	tights petite		
			right	left	right	left				
mediven plus	•		DT405	DT406	DT407	DT408			•	•
mediven elegance										
mediven for men										
mediven active										

mediven plus® PIP codes ccl 3

ccl 3 mediven plus PIP codes below knee									
		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT401/0/...	033-5984	323-6510	323-6528	323-6536	323-6544	323-6551	323-6569
	black	DT401/5/...	323-6643	323-6650	323-6668	323-6676	323-6684	323-6692	323-6700
open toe, petite	beige	DT402/0/...	323-6718	323-6726	323-6734	323-6742	323-6759	323-6767	323-6775
	black	DT402/5/...	323-6783	323-6791	323-6809	323-6817	323-6825	323-6833	323-6841
open toe, extra wide	beige	DT421/0/...	323-6577	323-6585	323-6593	323-6601	323-6619	323-6627	323-6635

ccl 3 mediven plus PIP codes thigh length

Silicone topband		medi code	I	II	III	IV	V	VI	VII
open toe	beige	DT400/0/...	323-9902	323-9910	323-9928	323-9936	323-9944	323-9951	323-9969
	black	DT400/5/...	323-9993	324-0009	324-0017	324-0025	324-0033	324-0041	324-0058
open toe, thigh extra wide	beige	DT400X/0/...	334-2078	334-2086	334-2094	334-2110	334-2128	334-2144	334-2151
open toe, petite	beige	DT399/0/...	323-9340	323-9357	323-9365	323-9373	323-9381	323-9399	323-9407
	black	DT399/5/...	323-9415	323-9423	323-9431	323-9449	323-9456	323-9464	323-9472
open toe, thigh extra wide, petite	beige	DT399X/0/...	334-2250	334-2268	334-2292	334-2300	334-2318	334-2326	334-2334

ccl 3 mediven plus PIP codes thigh length									
Waist attachment		medi code	I	II	III	IV	V	VI	VII
open toe, right	beige	DT405/0/...	324-0223	324-0231	324-0249	324-0256	324-0264	324-0272	324-0280
	black	DT405/5/...	324-0298	324-0306	324-0314	324-0322	324-0330	324-0348	324-0355
open toe, petite, right	beige	DT407/0/...	323-9621	323-9639	323-9647	323-9654	323-9662	323-9670	323-9688
	black	DT407/5/...	323-9696	323-9704	323-9712	323-9720	323-9738	323-9746	323-9753
open toe, left	beige	DT406/0/...	324-0082	324-0090	324-0108	324-0116	324-0124	324-0132	324-0140
	black	DT406/5/...	324-0157	324-0165	324-0173	324-0181	324-0199	324-0207	324-0215
open toe, petite, left	beige	DT408/0/...	323-9480	323-9498	323-9506	323-9514	323-9522	323-9530	323-9548
	black	DT408/5/...	323-9555	323-9563	323-9571	323-9589	323-9597	323-9605	323-9613

UCS™ Debridement

Sterile, pre-moistened debridement cloth.

- effective, gentle debridement
- removes wound debris, necrotic material, slough, exudate
- effective for removal of hyperkeratotic tissue
- cleans and moisturises the whole limb

UCS Debridement is a class IIb sterile device that is recommended for use with chronic and acute wounds, ulcers of all types, pressure sores and 1st or 2nd degree burns. It acts immediately, does not inhibit granulation and is compatible with subsequent use of any type of dressing.



UCS Debridement

UCS Debridement

medi code	DT500
PIP code	384-4271

For the latest articles by clinicians on UCS Debridement, email us at UCS@mediuk.co.uk

Wound debridement at an early and appropriate stage is likely to accelerate wound healing and improve patient care. Effective debridement has been shown to be associated with reduced exudate, reduced odour and the appearance of granulation in the wound bed.*

The UCS Debridement has three specific properties that combine for effective debridement:

- it is supplied pre-moistened with an anallergic solution to soften the necrotic tissue, plaque, fibrin and slough
 - effective debridement straight from a packet
 - the solution aids cleaning and hydration of the whole limb without the need of a bucket of water
- UCS Debridement is supplied in individual sterile foil packs in boxes of 10. Code no. DT500



Out of bandages



After one UCS in less that 4 minutes

*Downe A. (2012): How wound cleansing & debriding aids management and healing.

juxtacures®

Compression Ulcer Recovery System

The simple and effective alternative to bandages for treatment of venous leg ulcers.

Features

- ability to reduce size of garment as the limb volume decreases
- Breathe-O-Prene® fabric: breathable, flexible and comfortable material
- anti-odour properties
- 6 months sustained compression
- machine washable
- lightweight and discreet
- Accurate pressure measurement with Built-In-Pressure system (BPS) for safe application of compression

Benefits

- effective measurable and instantly adjustable compression
- effective compression reduces exudate levels
- cost effective
- saves time and reduction in clinical hours/clinical waste
- easy donning and doffing promotes patient empowerment
- allows patient to self-manage compression and hygiene needs
- allows free range of motion at the ankle to encourage calf muscle pump activity
- encourages the use of normal footwear

For your copy of the NICE Medtech Innovation Briefing,
email enquiries@mediuk.co.uk




Product benefits



CIRCAID®

Measurement guide



L: length of garment following contour of leg

Comfort Compression Anklet

Standard **Large**

Arch circumference (εH) 19 – 36 cm 22 – 44 cm

juxtacures length measurements

	Short	Standard	Long
L	28 cm	33 cm	38 cm
medi code	DT 100/2	DT 100/3	DT 100/4
PIP code	368-0733	368-0741	368-0758

Extras available:

Comfort liner (contains 2 per pack)

Comfort Compression Anklet (contains 2 per pack)
standard and large sizes

juxtacures extras

	Comfort liner	Comfort Compression Anklet Standard	Large
medi code	DT 100L	DT 100SA	DT 100LA
PIP code	377-1086	377-1102	377-1110

juxtalite[®] lower leg

Instantly adjustable inelastic compression device for venous diseases in patients who cannot tolerate or apply compression garments.

Features

- Accurate pressure measurement with Built-In-Pressure system (BPS) for safe application of compression
- inelastic material for high working pressure and low resting pressure
- garment straps can be individually readjusted while on the limb
- anti-odour and anti-microbial fabric
- latex free
- sustained compression with defined mmHg for up to 24 hours a day if required
- available in 8 sizes and 2 lengths



Product benefits




CIRCAID[®]

pack contains

- 1 juxtalite
- 1 BPS card
- 2 Comfort Liners
- 2 Comfort Compression Anklets

For the latest articles by clinicians on juxtalite, email us at enquiries@mediuk.co.uk

Measurement guide



Circumferences in cm

	S	M	M-XW	L
e C	26-36	33-44	44-54	41-51
e B	19-22	22-27	22-27	25-32
	L-XW	XL	XL-XW	XXL
e C	51-61	37-48	48-61	53-64
e B	25-32	30-37	30-37	35-42

Actual garment length in cm

	Short	Long
l B-D	28	33

All garments contain the circaid Built-In-Pressure system (BPS).

juxtalite[®] lower leg PIP codes

juxtalite, standard profile						
Length	medi code	small (2)	medium (3)	large (4)	x large (5)	xx large (6)
short 28 cm	DT301/...	394-0996	394-1002	394-1010	394-1028	394-1036
long 33 cm	DT300/...	394-0913	394-0921	394-0939	394-0947	394-0954

juxtalite, extra wide profile				
Length	medi code	medium xw (3)	large xw (4)	x large xw (5)
short 28cm, extra wide	DT301X/...	394-1044	394-1051	394-1069
long 33 cm, extra wide	DT300X...	394-0962	394-0970	394-0988

NICE Medtech Innovation Briefing supports the use of juxtacures for treating venous leg ulcers.

A recently published NICE Medtech Innovation Briefing (MIB25) in March 2015 suggests that juxtacures is an effective method of compressing venous leg ulcers and that it shows considerable cost savings compared to conventional bandaging techniques.

The juxtacures is simple to apply and adjust and the pressure applied can be easily checked by a clinician, carer or health care professional.

Effective compression with juxtacures leads to **reduced exudate** and therefore **less reliance on absorbent dressings**. The NICE MIB (Medtech Innovation Briefing) also highlights the **reduction in clinician's time** that can be realised when juxtacures is used, both through a reduced number of visits per week and a reduced duration

of visit due to the ease of removal and re-application. The document features a recent case study (Bianchi et al. 2013) which reported **cost savings with the juxtacures compared with conventional bandaging** over a 6 month period. It illustrated cost savings in these three areas:

- dressings: average saving of £753 per patient
- bandages: average saving of £881 per patient
- Clinician time: average saving of £3,172 per patient.

The total average saving was £4,806 per patient.

Copies of the NICE Medtech briefing are available to download at www.mediuk.co.uk/juxtacures or hardcopies are available on request from medi UK email enquiries@mediuk.co.uk

mediven ulcer kit[®]

Cost effective treatment for crural venous ulcers.

Features

- ensures rapid application of accurate compression
- less bulky than bandaging, improving comfort and mobility
- offers an increase in patient confidence, independence and concordance
- machine washable
- liner provides a constant 20 mmHg, over stocking adds a further 20 mmHg whilst patient is active
- supplied with two liner stockings, one to wear and one to wash



Product benefits



Inspection opening



Soft cuff




Colour coded instep/heel area aids correct fitting of garment

Measurement guide mediven ulcer kit®

mediven ulcer kit				
Size	Leg measurements (circ. in cm)			Colour mark
	Ankle ϵ B	Calf ϵ C	Below knee ϵ D	Heel
I	18 – 20	28 – 34	27 – 33	pink
II	20 – 22	30 – 37	29 – 36	blue
III	22 – 24	33 – 40	32 – 39	yellow
IV	24 – 26	35 – 43	34 – 42	red
V	26 – 28	37 – 46	36 – 45	grey
VI	28 – 30	39 – 49	38 – 48	green
VII	30 – 32	41 – 51	40 – 50	brown

Length: Floor (A) to below knee (D) in cm		
	Petite	Regular
A – D calf stockings	34 – 38	39 – 44
mediven ulcer kit	DT 276	DT 275
liner refill (2 per pack, universal length)	–	DT 277



mediven ulcer kit® PIP codes

Ordering information

The mediven ulcer kit can be ordered as a kit or as individual components.

The kit contains:

- 2 x liners (understockings)
- 1 x mediven ulcer plus over stocking
- 1 x easy-on applicator

Additional mediven plus stockings and under stockings can be ordered separately.

mediven ulcer kit PIP codes	medi code	I	II	III	IV	V	VI	VII
mediven ulcer kit, standard	DT275/...	329-2109	329-2117	329-2125	329-2133	329-2141	329-2158	329-2174
mediven ulcer kit, petite	DT276/...	329-2190	329-2208	329-2216	329-2224	329-2232	329-2240	329-2257
Under stocking, liner pack (2 x)	DT277/...	329-2265	329-2273	329-2281	329-2299	329-2307	329-2315	329-2323

mediven ulcer kit® – for new ulcers or small non exudating ulcers

Haemodynamic effects with guaranteed success

In addition to resting and working pressures, material properties such as stiffness and hysteresis are particularly important for effective compression therapy. Studies have confirmed the outstanding effects of the mediven ulcer kit in preventing and reducing oedema. The unique force/stretch properties and high dynamic stiffness values optimise exercise therapy for crural ulcers.

Plan the treatment results with the patient

mediven ulcer kit is simple to use, is very comfortable to wear, and is eminently suitable for use in everyday life. This results in high levels of acceptance

among patients who consequently accept personal responsibility for their own treatment. Two compressive components ensure efficient compression therapy and extend effective treatment time by up to eight hours a day.

Treat economically – preserve resources

mediven ulcer kit is clearly superior to traditional treatment approaches:

- quick and simple donning with little effort
- enables healthcare assistants to apply accurate graduated compression with the highest medical efficiency
- simple wound inspection
- shorter treatment time



Lymphoedema

What is oedema?

Oedema is a visible and palpable accumulation of fluid in the tissue. Experts distinguish between localised oedema (e.g. swelling of the eyelids, swelling of an arm or leg) and generalised oedema, in which the entire body is swollen. The most clinically important of the numerous forms of oedema are lymphoedema, lipoedema and phleboedema.

Lymphoedema

When the lymphatic system can no longer fulfil the function described above, lymphoedema develops. Lymphoedema consists of collections of fluid in the tissues, mainly in the limbs and rarely throughout the body.

Signs of lymphoedema are

- gradual or sudden swelling of the affected limb, rarely of a number of limbs
- taut skin and a feeling of heaviness at the affected site
- positive Stemmer's sign – the skin on the second toe cannot be raised with two fingers and feels taut and hard
- the skin pits when pressed with a finger at the affected site
- the natural skin folds on the toes and fingers become much deeper
- at later stages hardening (fibrosis) of the affected site

Diagnosis

Only an expert can determine

whether a patient actually has lymphoedema. Usually the specialist can recognise which type of oedema it is from the patient's history and by inspection and palpation (palpating the surface of the body).

If the clinical picture is unclear and in order to rule out malignant processes, modern medical imaging technology may be used in some cases.



Positive Stemmer's sign

Lymphoedema – Types

There are basically two types of lymphoedema.

Primary lymphoedema

is congenital. In this case, some areas of the lymphatic system show defective development. Although damage is present, there is often little or no sign of primary lymphoedema before puberty and it is only when there is a change in the hormone levels, causing the lymph-obligatory load to increase to such an extent that the limited capacity of the lymphatic system to transport it away is no longer adequate, that an oedema develops and eventually becomes evident.

Secondary lymphoedema

is the result of various conditions:

- after removal of lymph nodes during cancer surgery
- damage to the lymphatic system as a result of surgical interventions
- radiation
- infections with bacteria, fungi or tropical nematodes (filariasis)
- others (injuries, burns, bone fractures)
- lipoedema

Lymphoedema – Stages

Stage 0

As yet unidentified damage to the lymphatic vessels, no lymphoedema.

Stage 1 – reversible stage

The (soft) oedema appears over the course of the day and is reduced or disappears entirely when the limbs are elevated. The tissue pits when pressed and the area holds the indentation for a short time.



Stage 1

Stage 2 – spontaneously irreversible stage

The oedema persists even after prolonged periods of rest. The skin is often already fibrotic and elevation does not reduce swelling. There is little or no pitting when the skin is pressed.

Stage 3 – elephantiasis

Complicated swelling characterises by gross skin changes and fibrotic tissue, which may be wart-like or take the form of small blisters or fistulas, from which lymph leaks out, or can reach the point of huge swellings (“elephantiasis”).



Stage 2



Stage 3



Functional and fashionable

mediven® mondi flat knit compression available on the Drug Tariff.

Fresh new colours without compromising lymphoedema treatment.

Choose from aqua, medi magenta, navy, grey, black, sand, caramel, cashmere, moss-green and cherry-red.

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

mediven mondi® made-to-measure

Flat knit made-to-measure garment
for effective control of lymphoedema.

Features

- made-to-measure options on Drug Tariff: below knee, thigh length, thigh length with waist attachment, one legged panty and tights
- simple measuring form – available on request
- wide range of colours available (caramel, sand, navy, black, grey, magenta, aqua, cashmere, moss-green, cherry-red)
- closed toe, silicone topband and zip variants available in this range



Product benefits





Suitable for use with this product
and available on prescription:
medi 2in1

helps to put on and take off open toe
and closed toe compression garments.
Product code DT904, PIP code 329-8932



Measurement guide mediven mondi®



ℓ = length measurements
 c = circumferences

mediven mondi made-to-measure on FP10 is available in below knee, thigh length, thigh length with waist attachment, one legged panty and tights. Please contact medi customer services for your free measurement pad and guide.

RAL compression class 1

Suitable for treatment of mild lymphoedema, mild varicosis without oedema and for early varicosis in pregnancy.


Colours for compression class 1 leg garments / 18 – 21mmHg						
Styles	open toe	closed toe	AD	AG / SN	AG	ATE / ATEH
						
			calf	thigh with or without topband	thigh with waist attachment right left	one legged panty
mediven mondi made-to-measure	•	•	DTMK1 DTMT1		DTMW1	DTMOP1

Product codes for compression class 1 leg garments / 18 – 21mmHg

colours at no extra cost



mediven mondi
made-to-measure



mediven mondi® codes ccl 1

ccl 1 mediven mondi codes

	medi code
below knee	DTMK1
thigh length	DTMT1
thigh length with waist attachment	DTMW1
one legged panty	DTMOP1

Add

	medi code
silicone topband	DTTB
zip	DTZIP
closed toe	DTTOE

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

RAL compression class 2

Suitable for mild/moderate lymphatic conditions, pronounced varicosis with oedema present, after healing moderate ulcers.

Product codes for compression class 2 leg garments / 23 – 32mmHg								
Styles	open toe	closed toe	AD	AG / SN	AG	AT	ATE / ATEH	
								
			calf	thigh with or without topband	thigh with waist attachment right left	tights (pantyhose)	one legged panty	
mediven mondi made-to-measure	•	•	DTMK2	DTMT2	DTMW2	DTMTG2	DTMOP2	

colours at no extra cost

Product codes for compression class 2 leg garments / 23 – 32mmHg



mediven mondi
made-to-measure

-
-
-
-
-
-
-
-
-
-
-

mediven mondi® codes ccl 2

ccl 2 mediven mondi codes

	medi code
below knee	DTMK2
thigh length	DTMT2
thigh length with waist attachment	DTMW2
tights	DTMTG2
one legged panty	DTMOP2

Add

	medi code
silicone topband	DTTB
zip	DTZIP
closed toe	DTTOE

product information

measurement guide

round knit ccl 1

round knit ccl 2

round knit ccl 3

wounds

oedema

RAL compression class 3

Suitable for lymphoedema, significant varicosis, after healing long standing ulcers or if condition fails to respond to class 2 compression.

Product codes for compression class 3 leg garments / 34 – 46mmHg							
Styles	open toe	closed toe	AD	AG / SN	AG	AT	ATE / ATEH
							
			calf	thigh with or without topband	thigh with waist attachment right left	tights (pantyhose)	one legged panty
mediven mondi made-to-measure	•	•	DTMK3	DTMT3	DTMW3	DTMTG3	DTMOP3

Product codes for compression class 3 leg garments / 34 – 46mmHg

colours at no extra cost

beige black navy grey caramel sand aqua medi magenta cashmere moss-green cherry-red



mediven mondi
made-to-measure

• • • • • • • • • • •

mediven mondi® codes ccl 3

ccl 3 mediven mondi codes

	medi code
below knee	DTMK3
thigh length	DTMT3
thigh length with waist attachment	DTMW3
tights	DTMTG3
one legged panty	DTMOP3

Add

	medi code
silicone topband	DTTB
zip	DTZIP
closed toe	DTTOE

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema



How wearers can benefit from the comfort of flat knit

“I developed Lymphoedema at the age of 27, 4 years on I was still coming to terms with it.

It happened overnight when my right leg appeared to be swollen. I could not recall hurting it or spraining anything it just seemed to suddenly and painlessly appear. After many tests and scans and Google searches I was diagnosed with Lymphoedema.

I was prescribed stockings and taught Simple Lymphatic Drainage massage (SLD) which helps keep the oedema under control. I wear mediven mondi stockings, they are flat knit which I was unsure of at first but they look great with a short skirt which is something I never thought I would wear again.”

Donna Hughes

mediven esprit® made-to-measure

Flat knit made-to-measure garment
for effective control of lymphoedema.

Features

- made-to-measure options on Drug Tariff
- simple measuring form – available on request
- wide range of colours available (caramel, sand, navy, black, grey, magenta, aqua, cashmere, moss-green, cherry-red)
- shoulder cap, silicone topband, zip and pattern variants available in this range



Product benefits



Suitable for use with this product
and available on prescription:

medi 2in1

helps to put on and take off open toe
and closed toe compression
garments. Product code DT904, PIP
code 329-8932



product
information

measurement
guide

round knit
ccl 1




round knit
ccl 2

round knit
ccl 3

wounds

oedema

Measurement guide mediven esprit®



l = length measurements
 e = circumferences

mediven esprit made-to-measure on FP10 is available in armsleeve, combined armsleeve with handpiece, combined armsleeve with fingers, handpiece and glove. Optional extras are shoulder cap, silicone topband and zip. Please contact medi customer services for your free measurement pad and guide.

RAL compression class 1

Suitable for treatment of arm lymphoedema

Product codes for compression class 1 arm garments / 18 – 21mmHg

Styles

armsleeve



combined
armsleeve with
handpiece



combined
armsleeve with
fingers



handpiece



glove



mediven esprit
made-to-measure

DTECG1

DTEAG1

DTEAF1

DTEHP1

DTEGV1

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

RAL compression class 2

Suitable for treatment of arm lymphoedema

Product codes for compression class 2 arm garments / 23 – 32 mmHG

Styles

armsleeve



combined
armsleeve with
handpiece



combined
armsleeve with
fingers



handpiece



glove



mediven esprit
made-to-measure

DTECG2

DTEAG2

DTEAF2

DTEHP2

DTEGV2

RAL compression class 3

Suitable for treatment of arm lymphoedema

Product codes for compression class 3 arm garments / 34 – 46 mmHg

Styles

armsleeve



combined
armsleeve with
handpiece



combined
armsleeve with
fingers



handpiece



glove



mediven esprit
made-to-measure

DTECG3

DTEAG3

DTEAF3

DTEHP3

DTEGV3

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

colours at no extra cost

Product codes for compression class 1, 2 & 3 arm garments



mediven esprit
made-to-measure

-
-
-
-
-
-
-
-
-
-
-

mediven esprit® codes ccl 1, ccl 2 & ccl 3

	ccl 1 mediven esprit codes			medi code		
	ccl 1	ccl 2	ccl 3			
armsleeve	DTECG1	DTECG2	DTECG3			
combined armsleeve with handpiece	DTEAG1	DTEAG2	DTEAG3			
combined armsleeve with fingers	DTEAF1	DTEAF2	DTEAF3			
handpiece	DTEHP1	DTEHP2	DTEHP3			
glove	DTEGV1	DTEGV2	DTEGV3			

	Add	medi code
silicone topband		DTSTB
zip		DTZIP
shoulder cap		DTESC



Live, Laugh, Love



Sportive



Timeless

juxtafit® armsleeve and hand wrap

Instantly adjustable inelastic compression device for a cost effective alternative to bandaging for the control of arm lymphoedema.

Features

- inelastic adjustable gradient compression
- dual use for day and/or night
- comfort cover up sleeve included
- 6 standard sizes and 2 lengths – right and left side specific
- six month warranty
- high working pressure
- low resting pressure
- handwrap also available
- easy to self-manage



Product benefits



CIRCAID®

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

Measurement guide juxtafit® armsleeve

juxtafit armsleeve circumference measurements: size in cm


	1	2	3	4	5	6
c G top of upper arm	28–32	30–34	32–37	35–40	38–44	42–48
c F mid upper arm	26–30	28–32	30–35	33–38	36–41	39–45
c E elbow	25–28	27–30	29–32	31–34	33–36	36–40
c D mid forearm	23–26	25–28	27–30	29–32	31–34	33–37
c C wrist	15–17	17–19	19–21	21–23	23–25	25–28

juxtafit armsleeve length measurements: size in cm


	Standard	Long
l C – E wrist to elbow	up to 25	over 25

juxtafit hand wrap circumferences (cB): size in cm

XS	S	M	L	XL
14–17	17–20	20–23	23–26	26–29



juxtafit armsleeve



juxtafit hand wrap

juxtafit® armsleeve order and PIP codes

juxtafit armsleeve PIP codes								
		medi code	1	2	3	4	5	6
Armsleeve right	standard length	DT970R/...	366-9132	366-9140	366-9157	366-9165	366-9173	366-9181
	long	DT980R/...	366-9256	366-9264	366-9272	366-9280	366-9298	366-9306
Armsleeve left	standard length	DT970L/...	366-9199	366-9207	366-9215	366-9223	366-9231	366-9249
	long	DT980L/...	366-9314	366-9322	366-9330	366-9348	366-9355	366-9363

juxtafit hand wrap PIP codes							
		medi code	1	2	3	4	5
Hand wrap right		DT971R/...	366-9371	366-9389	366-9397	366-9405	366-9413
Hand wrap left		DT971L/...	366-9421	366-9439	366-9447	366-9454	366-9462

mediven harmony® & mediven harmony® short

The revolutionary round knit armsleeve.
Now in a short variant.

Features

- double ellipse around elbow
- 2 lengths
- extra wide sizes
- available in sizes 1 – 8 (I-VIII)

Benefits

- the double ellipses on the point and in the crease of the elbow offer comfort and freedom of movement
- prevents soreness and potential rubbing in elbow flexure
- stiffer fabric feels soft on the skin for comfort and clinical effectiveness




Product benefits



Double ellipse at the elbow for improved wearing comfort



sand available in regular armsleeve and glove



caramel available in short armsleeve

Measurement guide mediven harmony®

Variants (all variants available as regular or extra wide top)

AG combined armsleeve



Regular CCL1: DT734H
CCL2: DT736H
Extrawide CCL1: DT735H
CCL2: DT737H

AG combined armsleeve with grip top



Regular CCL1: DT754H
CCL2: DT758H
Extrawide CCL1: DT755H
CCL2: DT759H

AG combined armsleeve with shouldercap



Regular CCL1: 740H
CCL2: 744H
Extrawide CCL1: 741H
CCL2: 745H

CG



Regular CCL1: DT714H
CCL2: DT716H
Extrawide CCL1: DT715H
CCL2: DT717H

CG with grip top



Regular CCL1: DT746H
CCL2: DT750H
Extrawide CCL1: DT747H
CCL2: DT751H

CG with shouldercap



Regular CCL1: 738H
CCL2: 742H
Extrawide CCL1: 739H
CCL2: 743H

Variants

AC₁ glove with compressive fingers



Regular CCL1: DT760H
CCL2: DT761H

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

Measurement guide mediven harmony® short

AG combined armsleeve with grip top



Regular	CCL1: DT754HS CCL2: DT758HS
Extrawide	CCL1: DT755HS CCL2: DT759HS

CG with grip top



Regular	CCL1: DT746HS CCL2: DT750HS
Extrawide	CCL1: DT747HS CCL2: DT751HS


Variants

AC₁ glove with compressive fingers



Regular	CCL1: DT760H CCL2: DT761H
----------------	------------------------------

Measurement guide mediven harmony® glove



Length (cm) all sizes



AB	3.5 – 4.5
AC	10 – 12
AC ₁	15 – 16

Circumference / glove (cm)

Sizes	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
cC1 = 100 % *	15 – 18	17 – 20	19 – 21.5	21 – 23.5	23 – 26	25.5 – 30
cC	13.5 – 15.5	15 – 18	17 – 20	19 – 21	20 – 23.5	22.5 – 27.5
cB	15.5 – 20	17.5 – 21	19.5 – 23	20.5 – 24	22.5 – 26	24.5 – 30
cA	15 – 19	17 – 20	19 – 22	20 – 23	22 – 25	24 – 29

* Measurement cC1 is 5 cm above cC

Measurement guide mediven harmony® & mediven harmony® short armsleeves



Length (cm)

Standard length 41 – 47

Short length < 41

Circ (cms)	Size 1 (I)		Size 2 (II)		Size 3 (III)		Size 4 (IV)	
	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide
eG	25.5 – 28	28.5 – 31	28 – 30.5	30.5 – 33.5	30 – 33	33 – 36	32 – 35	35 – 38.5
eF	23.5 – 25.5	25.5 – 28	25 – 27.5	27.5 – 30	27.5 – 30	30 – 33	29.5 – 32	32 – 35
eE	23 – 25	24 – 26	24.5 – 26.5	26.5 – 29	26 – 28.5	28.5 – 31	27.5 – 30	30 – 33
eD	22 – 24	24 – 26	23 – 25	25 – 27.5	24 – 26	26 – 28.5	25 – 27.5	28 – 30.5
eC 1 = 100%*	15 – 16.5		16.5 – 18		18 – 19.5		19.5 – 21	
eC	13.5 – 14.5		14.5 – 15.5		16 – 17.5		17.5 – 19	
eA	15 – 18		16 – 19		17 – 20		19 – 22	
Circ (cms)	Size 5 (V)		Size 6 (VI)		Size 7 (VII)		Size 8 (VIII)	
	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide
eG	33.5 – 36.5	37 – 40.5	36 – 39.5	39.5 – 43	38 – 41.5	41.5 – 45.5	40.5 – 44.5	45 – 49.5
eF	31 – 34	34 – 37	33 – 36	36 – 39.5	35 – 38.5	38.5 – 42	37.5 – 41	41.5 – 45.5
eE	29 – 31.5	31.5 – 34.5	30.5 – 33.5	33.5 – 36.5	32.5 – 35.5	35.5 – 39	35 – 38.5	38.5 – 42
eD	27 – 29.5	29.5 – 32	28.5 – 31	31 – 34	30.5 – 33.5	33.5 – 36.5	33 – 36	36 – 39.5
eC 1 = 100%*	21 – 23		23 – 25		25 – 27.5		27.5 – 30	
eC	19 – 20.5		20.5 – 22.5		22.5 – 24.5		25 – 27.5	
eA	20 – 23		22 – 25		24 – 27		26 – 29	

*Measurement eC1 is 5 cm above eC

product information

measurement guide

round knit ccl 1

round knit ccl 2

round knit ccl 3

wounds

oedema

mediven harmony® ccl 1

mediven harmony ccl 1										
armsleeve	medi code		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	sand DT714H/...		381-7277	381-7285	381-7301	381-7319	381-7327	381-7335	381-7343	381-7368
CG xw	sand DT715H/...		381-7384	381-7392	381-7400	381-7418	381-7426	381-7434	381-7442	381-7459
CG with grip top	sand DT746H/...		381-7525	381-7533	381-7541	381-7558	381-7566	381-7574	381-7582	381-7590
CG with grip top xw	sand DT747H/...		381-7665	381-7673	381-7681	381-7699	381-7707	381-7715	381-7723	381-7731
combined armsleeve	medi code		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	sand DT734H/...		381-8002	381-8010	381-8028	381-8036	381-8044	381-8051	381-8069	381-8077
AG xw	sand DT735H/...		381-8085	381-8093	381-8101	381-8119	381-8127	381-8135	381-8143	381-8150
AG with grip top	sand DT754H/...		381-8168	381-8176	381-8184	381-8192	381-8200	381-8218	381-8226	381-8234
AG with grip top xw	sand DT755H/...		381-8242	381-8259	381-8267	381-8275	381-8283	381-8291	381-8309	381-8317
glove	medi code		2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)		
	sand DT760H/...		381-8549	381-8556	381-8564	381-8572	381-8580	381-8598		

mediven harmony® ccl 2

mediven harmony ccl 2										
armsleeve	medi code		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	sand DT716H/...		381-8606	381-8614	381-8622	381-8630	381-8648	381-8655	381-8663	381-8671
CG xw	sand DT717H/...		381-8689	381-8697	381-8705	381-8713	381-8721	381-8739	381-8747	381-8754
CG with grip top	sand DT750H/...		381-8762	381-8770	381-8788	381-8796	381-8804	381-8812	381-8820	381-8838
CG with grip top xw	sand DT751H/...		381-8846	381-8853	381-8861	381-8879	381-8887	381-8895	381-8903	381-8911
combined armsleeve	medi code		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	sand DT736H/...		381-9083	381-9091	381-9109	381-9117	381-9125	381-9133	381-9141	381-9158
AG xw	sand DT737H/...		381-9166	381-9174	381-9182	381-9190	381-9208	381-9216	381-9224	381-9232
AG with grip top	sand DT758H/...		381-9240	381-9257	381-9265	381-9273	381-9281	381-9299	381-9307	381-9315
AG with grip top xw	sand DT759H/...		381-9323	381-9331	381-9349	381-9356	381-9364	381-9372	381-9380	381-9398
glove	medi code		2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)		
	sand DT761H/...		381-9620	381-9638	381-9646	381-9653	381-9661	381-9679		

product information

measurement guide

round knit ccl 1

round knit ccl 2

round knit ccl 3

wounds

oedema

mediven harmony® ccl 1

armsleeve with shoulder cap		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	sand	DT738H/...	381-7830	381-7848	381-7855	381-7863	381-7871	381-7889	381-7897	381-7905
CG xw	sand	DT739H/...	381-7921	381-7939	381-7947	381-7954	381-7962	381-7970	381-7988	381-7996
combined armsleeve with shoulder cap		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	sand	DT740H/...	381-8325	381-8333	381-8341	381-8358	381-8366	381-8374	381-8382	381-8390
AG xw	sand	DT741H/...	381-8408	381-8416	381-8424	381-8432	381-8440	381-8457	381-8465	381-8473

mediven harmony® ccl 2

armsleeve with shoulder cap		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	sand	DT742H/...	381-8929	381-8937	381-8945	381-8952	381-8960	381-8978	381-8986	381-8994
CG xw	sand	DT743H/...	381-9000	381-9018	381-9026	381-9034	381-9042	381-9059	381-9067	381-9075
combined armsleeve with shoulder cap		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	sand	DT744H/...	381-9406	381-9414	381-9422	381-9430	381-9448	381-9455	381-9463	381-9471
AG xw	sand	DT745H/...	381-9489	381-9497	381-9505	381-9513	381-9521	381-9539	381-9547	381-9554

mediven harmony® short ccl 1

armsleeve with grip top		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	caramel	DT746HS/...	400-6631	400-6656	400-6664	400-6672	400-6680	400-6698	400-6706	400-6714
CG xw	caramel	DT747HS/...	400-6649	400-6722	400-6730	400-6748	400-6755	400-6763	400-6771	400-6789
combined armsleeve with grip top		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	caramel	DT754HS/...	400-6953	400-6979	400-6987	400-6995	400-7001	400-7019	400-7027	400-7035
AG xw	caramel	DT755HS/...	400-6961	400-7043	400-7050	400-7068	400-7076	400-7084	400-7092	400-7100

mediven harmony® short ccl 2

armsleeve with grip top		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
CG	caramel	DT750HS/...	400-6797	400-6813	400-6821	400-6839	400-6847	400-6854	400-6862	400-6870
CG xw	caramel	DT751HS/...	400-6805	400-6888	400-6896	400-6904	400-6912	400-6920	400-6938	400-6946
combined armsleeve with grip top		medi code	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)	8 (VIII)
AG	caramel	DT758HS/...	400-7118	400-7134	400-7142	400-7159	400-7167	400-7175	400-7183	400-7191
AG xw	caramel	DT759HS/...	400-7126	400-7209	400-7217	400-7225	400-7233	400-7241	400-7258	400-7266


The Whitaker pouch from medi

For treatment of scrotal oedema in palliative or acute situations.

- available in 3 sizes of pouch
- 100 cm belt as standard
- 2 lengths of belt extension available

Checking and refitting

- you are advised to teach your patient or their carers to check the skin condition regularly
- due to the possibility of rapid reduction of the oedema with this garment you may need to re-measure for a smaller garment more frequently than standard garments
- regular reviews are recommended



Fitting: The pouch is designed to fit the scrotum only, but may include the penis for palliative patients.

Product benefits



Measurement guide Whitaker pouch

Whitaker Pouch			
Circumference of scrotum	Size	medi code	PIP code
33 – 40 cm	Small	DT700/2	363-9713
41 – 48 cm	Medium	DT700/3	363-9721
49 – 56 cm	Large	DT700/4	363-9739

Waist circumference	Standard belt 1025402	40 cm extension required 1025404	80 cm extension required 1025401
100 – 150 cm	•		
145 – 190 cm	•	•	
185 – 225 cm	•		•

juxtafit® essentials upper leg with knee piece

Instantly adjustable inelastic compression device for effective control of oedema in the knee and thigh.

Features

- graduated profile
- alternating band system for easy application
- Breathe-O-Prene® for a flexible, comfortable fit
- low resting pressure & high working pressure
- compatible with juxtafit lower leg garments
- accurate pressure measurement with Built-In-Pressure system (BPS) for safe application of compression
- available in 5 sizes and 2 lengths. Left and right specific
- six month warranty
- full leg liner
- makes self-management easy



Product benefits



CIRCAID


juxtafit® essentials upper leg with knee order and PIP codes

Product Order Codes & PIP Codes:							
Sizes		medi code	XS (1)	S (2)	M (3)	L (4)	XL (5)
Short length (45 cm)	left	DT 955L/...	378-6498	378-6506	378-6514	378-6522	378-6530
	right	DT 955R/...	378-6548	378-6555	378-6563	378-6571	378-6589
Long length (55 cm)	left	DT 956L/...	378-6597	378-6605	378-6613	378-6621	378-6639
	right	DT 956R/...	378-6654	378-6662	378-6670	378-6688	378-6696

juxtafit essentials extras	
Full leg liner	
code number	DT957
PIP code	382-4273

Measurement guide juxtafit[®] essentials upper leg with knee

juxtafit essentials upper leg with knee piece					
Circumference (in cm)	XS (1)	S (2)	M (3)	L (4)	XL (5)
e G	42 – 62	54 – 70	60 – 80	69 – 90	77 – 101
e E1 (5 cm above E)	37 – 48	42 – 55	50 – 65	58 – 76	66 – 87
e D	27 – 35	31 – 41	35 – 46	40 – 52	44 – 58
Garment length (in cm)					
	Short	Long			
ℓ D-G	45	55			



juxtafit® lower leg

Simple and effective instantly adjustable inelastic compression device for oedema control.

Features

- easy for self-management
- graduated profile
- alternating band system for easy application
- minimal adjustments
- Breathe-O-Prene® for a flexible, comfortable fit
- low resting pressure
- high working pressure
- XW calf variants available



Product benefits



CIRCAID®

accurate pressure measurement with Built-In-Pressure system (BPS) for safe application of compression

ankle foot pieces

For patients with resistant foot/ankle oedema.

Features

- simple effective inelastic containment garment
- ideal for self-care
- the cost effective solution for oedema management



juxtafit interlocking ankle foot wrap (open heel)



juxtafit ankle foot wrap (closed heel)

product information

measurement guide

round knit ccl 1

round knit ccl 2

round knit ccl 3

wounds

oedema

Measurement guide juxtafit[®] lower leg & ankle foot wrap

juxtafit lower leg circumference: size in cm

Sizes	Small	Medium	Medium x-wide	Large	Large x-wide	X-Large	X-Large x-wide	XX-Large
C	30–38	36–43	43–50	41–48	48–56	41–51	51–61	53–64
B1	20–28	25–33	33–40	30–38	38–46	33–41	41–51	43–53
B	15–23	20–28	20–28	25–33	25–33	30–38	30–38	33–41


Product length: petite 28 cm or standard 36 cm (from B-D).

juxtafit ankle foot wrap available in 3 sizes and 2 styles sizes in cm


Sizes	Small	Medium	Large
H Arch circumference	19–24	24–29	29–34
B Ankle circumference around malleoli	<36	<36	<41

ankle foot wrap – **Closed heel**, foot length from base of toes > 16 cm

Interlocking ankle foot wrap – **Open heel**, foot length from base of toes 13–25 cm



juxtafit lower leg



juxtafit ankle foot wrap

juxtafit® lower leg & ankle foot wrap PIP codes

juxtafit lower leg PIP codes									
	medi code	Small	Medium	Medium x-wide	Large	Large x-wide	X-Large	X-Large x-wide	XX-Large
		2	3	3X	4	4X	5	5X	6
Petite length 28 cm (11")	DT950/...	348-0548	348-0555	368-0782	348-0563	368-0790	348-0571	368-0808	348-0589
Standard length 36 cm (14")	DT951/...	348-0597	348-0605	368-0816	348-0613	368-0824	348-0621	368-0832	348-0639

juxtafit ankle foot wrap PIP codes				
	medi code	Small	Medium	Large
		2	3	4
Open heel (interlocking)	DT960/...	348-0647	348-0654	348-0662
Closed heel	DT961/...	348-0670	348-0688	348-0696

duomed[®] soft

A British Standard (BS), soft elastic stocking listed as 'elastic hosiery' in Part IXA of the Drug Tariff for mild venous disease or palliative treatment.

Features

- sheer and soft fabric
- easy to apply
- 3 British Standard compression classes – ccl 1, ccl 2 & ccl 3 (black available in ccl 1 & ccl 2)
- below knee and thigh length available
- open or closed toe – compression class 1 & 2 (open toe only for class 3, in sand only)
- machine washable
- simple size selection
- colour: sand & black

Now available in black
(ccl 1 & ccl 2)

UNIQUE FEATURE: British Standard elastic hosiery with a silicone top band on the thigh length variant.



Thigh length with silicone top band

Circumference (in cm)					
Sizes	S	M	L	XL	XXL
⌒G Thigh	42–57	48–64	54–71	60–78	66–85
⌒C Calf	28–34	32–38	36–42	40–46	44–50
⌒B Ankle	19–21	22–24	25–27	28–30	31–34



Below knee

product
information

measurement
guide

round knit
ccl 1

round knit
ccl 2

round knit
ccl 3

wounds

oedema

duomed soft

duomed® soft PIP codes – sand colour

PIP codes for duomed soft – sand						
Sizes	medi code	S	M	L	XL	XXL
ccl 1						
below knee, open toe	DT511	367-4850	367-4868	367-4876	367-4884	367-4892
below knee, closed toe	DT512	367-4926	367-4934	367-4959	367-4967	367-4975
thigh length with silicone topband, open toe	DT611	367-4983	367-4991	367-5006	367-5014	367-5022
thigh length with silicone topband, closed toe	DT612	367-5030	367-5048	367-5055	367-5063	367-5071
ccl 2						
below knee, open toe	DT521	367-5089	367-5097	367-5105	367-5113	367-5121
below knee, closed toe	DT522	367-5139	367-5147	367-5154	367-5162	367-5170
thigh length with silicone top band, open toe	DT621	367-5188	367-5196	367-5204	367-5212	367-5220
thigh length with silicone topband, closed toe	DT622	367-5238	367-5246	367-5253	367-5261	367-5279
ccl 3						
below knee, open toe	DT532	367-5287	367-5295	367-5303	367-5311	367-5329
thigh length with silicone topband, open toe	DT632	367-5337	367-5345	367-5352	367-5360	367-5378

duomed® soft PIP codes – black colour

PIP codes for duomed soft - black						
Sizes	medi code	S	M	L	XL	XXL
ccl 1						
below knee, open toe	DT511/B/	388-1802	388-1810	388-1828	388-1836	388-1844
below knee, closed toe	DT512/B/	388-1851	388-1869	388-1877	388-1885	388-1893
thigh length with silicone topband, open toe	DT611/B/	388-1901	388-1919	388-1927	388-1935	388-1943
thigh length with silicone topband, closed toe	DT612/B/	388-1950	388-1968	388-1976	388-1984	388-1992
ccl 2						
below knee, open toe	DT521/B/	388-2008	388-2016	388-2024	388-2032	388-2040
below knee, closed toe	DT522/B/	388-2057	388-2065	388-2073	388-2081	388-2099
thigh length with silicone top band, open toe	DT621/B/	388-2107	388-2115	388-2123	388-2131	388-2149
thigh length with silicone topband, closed toe	DT622/B/	388-2156	388-2164	388-2172	388-2180	388-2198

Headquarters:

medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany

T +49 921 912-1371
F +49 921 912-8371
doctors@medi.de
www.medi-doctors.com

medi UK Ltd.
Plough Lane
Hereford HR4 OEL
Great Britain

T 01432 373500
F 01432 373510
enquiries@mediuk.co.uk
www.mediuk.co.uk

